



CrimpLoc® System



ILC Dover pioneered flexible containment for the pharmaceutical industry and continues to innovate processes for better performance. We invented the original CrimpLoc™ System to provide high containment and a secure closure when using the ILC Dover Continuous Liner, and other ArmorFlex® film solutions.

- One Tool for all Crimp Sizes
- Validated containment performance
- Light weight

Key Factors to Consider for Closure of in Process Pharmaceutical Packaging

- Containment
- Quality
- Ergonomics / Safety
- No plastic shards from separation cut
- Reliability



The CrimpLoc™ System is a complete kit to provide the secure closure for ArmorFlex® film packages. Kits are available including the CrimpLoc™ tool, Cutting Tool, and Crimps.

One CrimpLoc™ tool is needed for the range of ArmorFlex® closures offered by ILC Dover:

ArmorFlex® size	Crimp Designation
23" Diameter	Blue
14" Diameter	Red
4" Diameter	Yellow



CrimpLoc™ Kit



Crimp Removal Tool - The latest innovation for safe removal of the crimp. The crimp remains in one piece so no risk of contamination upon opening.

Beyond Boundaries™.

Innovators at our core, we develop engineered solutions for our customers' complex problems. Recognized globally for our flexible containment solutions, ILC Dover serves customers in a diverse range of industries, including pharmaceutical and biopharmaceutical manufacturing, personal care, food and beverage, chemical, aerospace, healthcare and government agencies. At ILC Dover, quality is a culture, not a measurement. Our customers will tell you that we cater to their every need and that we're highly innovative, responsive, dedicated and competitive. We have been innovating since 1947. ILC Dover's visionary solutions improve efficiency, safeguard workers and product, and prevent disasters – proof that we are on the front line of business excellence.

Engineering evolution beyond boundaries.



One Moonwalker Road
Frederica, DE 19946 USA
+1.302.335.3911
+1.800.631.9567

customer_service@ilcdover.com

www.ilcdover.com

